







## Solvent and chemicals resistant solution Format & standards

Solvents and chemicals resistant labels are recommended in histological staining. Labels can be immersed in different alcohol mixtures during dehydration, infiltration and tissue deposition in paraffin. They are resistant to the process of removing xylene and toluene from tissue fragments as well as staining and counter-staining with eosin, hematoxylin and other dyes. Label formats are matched to microscope slides and tissue cassettes which come in contact with solvents and chemicals. The labels do not need to be coated with additional transparent tape. They can also be stored for many years after being processed in archives.

SHAPE	SIZE (mm)	PRODUCTS
	23 x 19	
	23 x 23	

### Solution



Solvent and chemicals resistant solution

### These solutions are found in the following products



Labels for microscope slides